



TITLE: iSTAT1 Analyzer configuration Job Aid	Doc #: 13961
Section: \\Management System\PLM\Point of Care Testing\Multi Test Analyzer\	Version: 2.0 Current
Document Owner: Trifacility Supervisor	Effective Date: 7/18/2013
Final Approval: Dr Bassam Nassar	

1. Analyzer

Option	Description	Choices i-STAT1	Default0	ESMH
Page 1				
1- Language	Select a language for text	English, Japanese, Deutsch, Italiano, Nederlands, Espanol, Francais, Svensk or Portugueses	English	English
2- Date Format	For Clock Set function	mm/dd/yy dd/mm/yy	mm/dd/yy	dd/mm/yy
3- Sound	Optional beep when results or Quality Check are displayed	Disabled or Enabled	Enabled	Enabled
4- Auto-transmit	Analyzer automatically transmits results when placed in a Downloader	Disabled or Enabled	Enabled	Disabled
5- Memory Full	What will happen when analyzer's memory contains 5000 unsent results	Allow Test, Warn User or Lock Out	Allow Test	Allow Test
Page 2				
1-Batch Mode Timeout	Not active at this time			
2- Inactivity Timeout	Number of seconds after results are displayed and there is no operator intervention that analyzer will turn off	45 to 1620 seconds	<u>120 seconds</u>	240 sec
3- Upload Schedule	Time interval between uploads from analyzer to CDS Schedule Options if time interval exceeded	Disabled or Enabled with a specified time interval in hours from 1 to 65535 Allow Test, Warn User or Lock Out	Disabled	<input type="checkbox"/> Disabled
4- PCx Glucose Strip	Enables the PCx glucose test strip reader	Disabled or Enabled	Disabled	Disabled



Option	Description	Choices i-STAT1	Default0	ESMH
5-Clock Password	Password protection for the Set Clock function	Disabled or Enabled	Enabled	Enabled
Page 3				
1-Sync Clock	This will synchronize or update the real time clock in the analyzer to the CDS's clock at the time or each download. This option will eliminate the need to reset the analyzer's clock at the the beginning and end of Daylight Savings time.	Disabled or Enabled		Disabled

2. ID Entry (Operator ID)

Option	Description	Choices i-STAT1	Default0	ESMH
1- Operator ID				
1- Minimum Length	Minimum allowed ID length	Minimum = 0	0 digits	0 digits
2- Maximum Length	Maximum allowed ID length	Maximum = 15	<u>Maximum = 15</u>	<u>Maximum = 15</u>
3- Repeat ID	Analyzer can prompt for repeat entry of operator ID. This option can be enhanced by requiring the operator to repeat a manual or scanned entry	Disabled or Enabled Repeat Manual Entry Repeat Scanned Entry	Enabled	Disabled
4- Manual Entry	A check digit algorithm can be applied to manually entered operator and /or patient IDs	No Check Digit, Modulus 11 Check Digit or Modulus 10 Check Digit	No Check Digit	No Check Digit
5- Code I2 of 5	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled No Check Digit, USS or OPCC	<u>Enabled USS</u>	<u>Enabled USS</u>
Page 2				
1- Code 128	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	<u>Enabled</u>	<u>Enabled</u>
2- EAN-8,	Type of a valid barcode	Disabled or Enabled	<u>Enabled</u>	<u>Enabled</u>



Option	Description	Choices i-STAT1	Default0	ESMH
EAN-13	format for both the operator and <u>patient ID numbers</u>			
3- Codabar	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Enabled
4- Code 93	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Enabled
5- Code 39	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled Alphanumeric or Full ASCII	Enabled Full ASCII	Enabled Full ASCII
Page 3				
1- Code 39 Check Digit	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Disabled
2- Truncate First	User can select how to truncate digits from a scanned operator and/or patient ID	Enter number of leading characters to be stripped from the barcode.	0 Digits	0 Digits
3- Truncate Last	User can select how to truncate digits from a scanned operator and/or patient ID	Enter number of trailing characters to be stripped from the barcode.	0 Digits	0 Digits
4-Print ID	This option enables/disables printing of operator IDs on printouts from the Martel Printer. Disabling the printing of operator IDs can prevent uncertified operators from learning the IDs of certified operators.			Enabled



2. ID Entry (Patient ID)

Option	Description	Choices i-STAT1	Default0	ESMH
2- Patient ID				
1- Minimum Length	Minimum allowed ID length	Minimum = 0	0 digits	0 digits
2- Maximum Length	Maximum allowed ID length	Maximum = 15	Maximum = 15	Maximum = 15
3- Repeat ID	Analyzer can prompt for repeat entry of operator ID	Disabled or Enabled	Enabled	Disabled
4- ID Recall	Operator can recall last patient ID when analyzer prompts for patient ID	Disabled or Enabled	Enabled	Enabled
5- Manual Entry	A check digit algorithm can be applied to manually entered operator and /or patient IDs	No Check Digit, Modulus 11 Check Digit or Modulus 10 Check Digit	No Check Digit	No Check Digit
Page 2				
1- Code I2 of 5	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled No Check Digit, USS or OPCC	Enabled USS	Enabled USS
2- Code 128	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Enabled
3- EAN-8, EAN-13	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Enabled
4- Codabar	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Enabled
5- Code 93	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled	Enabled	Enabled
Page 3				
1- Code 39	Type of a valid barcode format for both the operator and patient ID numbers	Disabled or Enabled Alphanumeric or Full ASCII	Enabled Full ASCII	Enabled Full ASCII
2- Code 39	Type of a valid barcode	Disabled or Enabled	Enabled	Disabled



TITLE: iSTAT1 Analyzer configuration Job Aid	Version: 2.0 Current
---	-----------------------------

Option	Description	Choices i-STAT1	Default0	ESMH
Check Digit	format for both the operator and patient ID numbers			
3- Truncate First	User can select how to truncate digits from a scanned operator and/or patient ID	Enter number of leading characters to be stripped from the barcode.	0 Digits	0 Digits
4- Truncate Last	User can select how to truncate digits from a scanned operator and/or patient ID	Enter number of trailing characters to be stripped from the barcode.	0 Digits	0 Digits

3. Patient Tests

Option	Description	Choices i-STAT1	Default0	ESMH
1- Cartridge Auto-Chart	If enabled, the Chart page is displayed automatically when testing a cartridge and/or strip.	Disabled Enabled	Disabled	Disabled
2- Cartridge Information	If “Required information before running cartridge” is required, the operator must enter the operator and patient IDs before the analyzer will initiate a cartridge test cycle.	Not Required Required	Not_Required	Required
3-Cartridge Barcode		Not Required Required		Required
4- Cartridge Lot Number	If “Enter Lot Number” is required, the cartridge lot number must be entered before the analyzer will initiate a cartridge test cycle.	Not Required Required	Not Required	Required
5- Comment Code In Range	If required, the analyzer will prompt for a 1 to 3 character comment code after results for patient samples are displayed.	Disabled Optional Required	Disabled	Disabled
5- Comment Code Out of Range	If required, the analyzer will prompt for a 1 to 3 character comment code after results for patient	Disabled Optional Required	Disabled	Disabled



TITLE: iSTAT1 Analyzer configuration Job Aid	Version: 2.0 Current
---	-----------------------------

Option	Description	Choices i-STAT1	Default0	ESMH
	samples are displayed.			
Result Output				Disabled
Downloader Lockout				Disabled

1. QC Tests

Option	Description	Choices i-STAT1	Default0	ESMH
1-External Simulator	Time interval when a prompt for operator to use the external simulator should appear.	1- Disabled 2- Interval of specified hours (1 to 65535 hours)	Disabled	Disabled <u>744</u> hours
2- Internal Simulator	Time interval when the internal Electronic Simulator test will be run.	1- Disabled 2- 8/24 Hours (every 8 hours for blood gases and hematocrit and every 24 hours for other tests.) 3- Interval of specified hours (1 to 65535 hours) 4- Patient Tests -interval of specified number of patient tests	Interval 24 hours	8/24 Hours
3- Internal Simulator Schedule Option	Schedule for control options. The behavior of the analyzer if the simulator test fails can also be specified. If the Schedule Option Lockout is selected, the analyzer will continue to perform the simulator test and will continue to display "FAIL" on subsequent cartridges until the test passes. If Lockout is not selected, the simulator test will not be initiated again until next scheduled time.	Allow Test or Lock Out	Lock Out	Lock Out



TITLE: iSTAT1 Analyzer configuration Job Aid	Version: 2.0 Current
---	-----------------------------

5. Results

Note: Reference Ranges for pH, PCO2, PO2, lactate, and associated calculated results are venous ranges.

Option	Description		Choices i-STAT1	Default0	ESMH
1- Units and Ranges	Analyte Enabled or Disabled. When enabled, select reporting units for results.		Each test can be disabled. If disabled, the test will not appear on the results page and test results will not be stored. Select Units for analytes		
Analyte	Disabled ✓	Reference Range	Disabled or Enabled Unit Set	All tests enabled	
1- Cl		98 / 109	mmol/L or mEq/L	mmol/L	
2- K		3.5 / 4.9	mmol/L or mEq/L	mmol/L	mmol/L
3- Na		138 / 146	mmol/L or mEq/L	mmol/L	mmol/L
4- BUN UREA		8 / 26 17 / 56	1- BUN 2- UREA- mmol/L, mg/dL or g/L	BUN mg/dL	Disabled UREA- mmol/L
5- iCa		1.12 / 1.32	mmol/L or mEq/L or mg/dl	mmol/L	mmol/L
Page 2					
1- pH		7.310 / 7.410	Not expressed in units	N/A	Disabled
2- Glu		70 / 105	mmol, mg/dL or g/L	mg/dL	mmol
3- TCO2		24 / 29	mmol/L or mEq/L	mmol/L	mmol/L
4- AnGap		10 / 20	mmol/L or mEq/L	mmol/L	Disabled
5- PO2		80 / 105	kPa or mmHg	mmHg	Disabled
Page 3					
1- PCO2		41.0 / 51.0	kPa or mmHg	mmHg	Disabled
2- BE		-2 / 3	mmol/L or mEq/L	mmol/L	Disabled
3- HCO3		23 / 28	mmol/L or mEq/L	mmol/L	Disabled
4- sO2		95 / 98	% or None	%	Disabled
5- Hb		12 / 17	mmol/L, g/dL or g/L	g/dL	



Option	Description		Choices i-STAT1	Default0	ESMH
Page 4					
1- Hct		38 / 51	None or %PCV	%PCV	None
2- Lac		0.36 / 1.25	mmol/L, mg/dL or g/L	mmol/L	Disabled
3- Crea		0.6 / 1.3	mg/dL or umol/L	mg/dL	umol/L
4- ACT		79 / 149	Expressed only as WBT sec	WBT sec	Disabled
5-INR			Reference Value		Reference Value
6-cTnl			ng/ml		Disabled
2- Options					
1- Decimal Separator		Select comma (,) or period (.)	Comma or period	(.) Period	(.) Period
2- Test Selection		Requires the operator to select tests to be reported from a cartridge test panel.	Disabled or Enabled	Disabled	Disabled
3- Hematocrit		CPB options. Prompt or no prompt for CPB compensation when cartridge includes hematocrit sensor, and if no prompt, automatically Adjust or Do not adjust.	Prompt CPB, CPB Always or CPB Never K3 EDTA or K2EDTA/Heparin/None. (NaEDTA is included in this option and None means no anticoagulant.)	Prompt CPB K3 EDTA	<input type="checkbox"/> Prompt CPB <input type="checkbox"/> K3 EDTA <input type="checkbox"/> K2 EDTA <input type="checkbox"/> CPB Always <input type="checkbox"/> K3 EDTA <input type="checkbox"/> K2 EDTA



Option	Description		Choices i-STAT1	Default0	ESMH
		Reference anticoagulant used to calculate hematocrit result.			<input type="checkbox"/> CPB Never <input type="checkbox"/> K3 EDTA <input type="checkbox"/> K2 EDTA
Base Excess Calculation		Select Base Excess of Extracellular Fluid (BEecf) or Base Excess of Blood (BEb).	ecf or blood	ecf	ecf
ACT-C	Will allow results to run with Prewarm or Nonwarm mode	Prewarm (37) or Nonwarm (ambient temp)	Prewarm or Non-Prewarm	Prewarm	Disabled
ACT-K	Will allow results to run with Prewarm or Nonwarm mode	Prewarm(37) or Nonwarm (ambient temp)	Prewarm or Non-Prewarm	Prewarm	Disabled
Print Reference Ranges		Reference ranges can be printed with results	Disabled or Enabled	Disabled	Disabled

6. Password

Option	Description	Choices i-STAT1	Default0	ESMH
Password	For Clock Set, Change customization, and Utility functions protection	0-5 digits	No password	enter



7. Restore Factory Settings

Option	Description	Choices i-STAT1	Default0	ESMH
Restore Factory Settings	The customization profile can be reset to factory settings. DEFAULT0 indicates that the analyzer has factory settings. If the default is changed on the analyzer, the profile is listed as 00000000.	Yes or No	Analyzer is shipped with factory settings	no