



Clinitek Status QA Log

Final Approval: Dr Bassam Nassar

Month: _____

Year: _____

Serial Number: _____

Daily checks	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31
Clean sample table																															
Record quality issues as required																															
Initials																															

Monthly	
Clean white calibration bar	
Perform QC and record results	
Date and initials	

Comments/Issues:



Final Approval: Dr Bassam Nassar
